

Moves on the Sheet... Not on the Feet!

HEELIFT[®] Glide
SUSPENSION BOOT



Introducing the Heelift[®] Glide Suspension Boot

How will the Heelift[®] Glide save you time? It's simple. A boot that stays in place, moves with the patient, and is comfortable will eliminate the need for repeated adjustments and reapplications. A boot that effectively offloads and keeps the patient's heel safe from pressure means less time spent on costly treatments for pressure ulcers. And for those who have an existing pressure ulcer, the effective offloading and comfort mean you can concentrate on healing the wound. Heelift Glide represents a brand new chapter in the Heelift story.

Of course, the Heelift Glide still features the heel offloading that you've come to expect from a Heelift. Just like the original Heelift Suspension Boot, the Heelift Glide is ready to use right out of the bag or can be customised to meet your patient's needs.

We designed the Heelift Glide by incorporating years of feedback from clinicians and patients, making this the best Heelift yet. Find out for yourself.


Patents Pending

Because proper fit is crucial to the success of your prevention and treatment program, Heelift Glide is available in three sizes.



Petite	70 - 130 lbs	6 - 10 in	46 - 64 in
Standard	120 - 250 lbs	8 - 14 in	60 - 77 in
Bariatric	220 - 600 lbs	12 - 23 in	65 - 80 in

What makes the Heelift® Glide the right boot for your patient?



- Suspends the heel in space for a pressure-free environment.
- 40% lighter, cooler and less bulky than a pillow-style boot.
- Soft inner material cradles the calf, redistributing pressure over a wider area and providing interior grip.
- Fabric covering, tricot backing and innovative stitching ensures there are no rough, hard or sharp edges that can damage sensitive skin.

HEELIFT[®]
Glide
SUSPENSION BOOT

Simple. Proven. Solutions.



ORTHOTICS LTD
Phone: 01986 798120
Fax: 01986 798040
email: info@vmorthotics.co.uk
www.vmorthotics.co.uk